

# AVANTIK CSF<sup>20</sup>

## AUTOMATIC FILM COVERSAMPLER

The Avantik CSF20 Automatic Film Coversampler is designed to provide high quality, coverslipped slides in an easy-to-use, automated solution.



1. Xylene Pump System
2. Starter Basket Holder for Input Slides
3. Output Tray for Empty Slide Basket
4. Waste Bottle
5. Output Tray for Finished Slides Basket
6. End Basket Loading Station
7. Coverslipping Film Holder
8. User-Friendly LCD Display

### Product Features

- ❖ Coverslips up to 1,200 slides per hour
- ❖ No glues or mounting media needed
- ❖ Ideal for H&E, IHC and cytology slide coverslipping
- ❖ Large, clear front provides protection against contamination
- ❖ Displays process status and remaining capacity for slides, xylene and coverfilm

### User Benefits

- ❖ Integrated carbon filter for standalone use
- ❖ Intuitive menu navigation via touch display
- ❖ LED interior lighting for great interior visibility
- ❖ Rapid slide drying for immediate storage
- ❖ All components (such as coverslipping film) can be easily inserted and exchanged
- ❖ Minimal maintenance required, compared to traditional glass coverslippers
- ❖ Convenient windows in the doors allow for quick exchange of new slide baskets without interrupting coverslipping process

Peace of Mind,  
Guaranteed Every Day.  
Avantik, **A Higher Standard.**



<b>Part Number</b>	CSF20
<b>Slide Basket Size</b>	20 pcs
<b>Coverslipping Time</b>	~3 seconds per slide
<b>Coverslipping Medium Capacity</b>	Bottle of 500 ml
<b>Coverslipping Film (width x length)</b>	24 mm x 70 m
<b>Coverslipping Length</b>	40 - 60 mm
<b>Waste Bottle Capacity</b>	100 ml
<b>Dimensions (width x depth x height)</b>	730 x 310 x 545 mm
<b>Weight</b>	143 lbs
<b>Power Supply</b>	115 - 230 V / 50-60 Hz

Accessories	Quantity	Part Number
Tissue-Tek Film Coverslipping, 60 m	5/Case	SL6070
Pathotape Film Coverslipping, 70 m	1/Each	SL6400
Medite Film Coverslipping, 66 m	1/Each	SL6543
Cutting Blade	5/Box	SE3153
Charcoal Filter, Disposable	1/Each	SY6319
Basket for 20 Slides, Plastic	1/Each	GL5074
Waste Bottle, 100 ml	1/Each	GL5075
Xylene Bottle, 500 ml	1/Each	GL5076
Medite Stainer Compatible Hanger	1/Each	EQ6156
Epredia Gemini Stainer Compatible Hanger	1/Each	EQ6157
Leica Stainer Compatible Hanger	1/Each	EQ6158



**SCAN NOW**

to learn more about  
the Avantik CSF20

Contact Us Today!

[avantik-us.com](http://avantik-us.com)

800.783.9424

