

# CX33 MICROSCOPE



The CX33 offers improved performance and increased cost-efficiency. It is ergonomically designed with low positioned coarse & fine focusing controls and convenient hand grips at the front and back. This microscope is intended for educational and medical applications.

## Product Information

**CX33** Olympus Microscope Systems

### Features

- LED illumination
- Trinocular observation tube
- 10X widefield eyepieces
- 2x, 4x, 10x, and 40x objectives
- 4 position nose-piece turret
- Infinity optical system

### Benefits

- Inward-facing quadruple revolving nosepiece
- Fixed abbe condenser
- Right hand stage
- 2.4w LED illumination

### Optional Accessories

- Dual-view side-by-side adapter
- Dual-view face-to-face adapter
- Digital camera mount
- LED monitor

Call us today.

800-783-9424  
info@avantik-us.com  
avantik-us.com

