

SAVE YOUR SPECIMEN BY ELIMINATING REAGENT MIX-UPS.

Match Before You Process.

Avantik Reagent Labels



ReOrder# R54029 32 Commerce Street, Springfield, NJ 07081
t 888.392.8411 www.avantik-us.com

ALCOHOL 100%



**Reagent
100% Alcohol**
200 PROOF REAGENT GRADE/ACS

Formula (Denatured 3A)	200 Proof	Percent by Volume
SDA 3A (Ethanol-Methanol)	95 parts by volume	Ethanol 90.48%
IPA	5 parts by volume	Methanol 4.52%
		IPA 5.00%

1 Gallon | Case = 4 x 1 Gallon

DANGER! Highly flammable liquid and vapor. Harmful if swallowed. Causes skin irritation. Causes serious eye irritation. May cause respiratory irritation. Causes damage to organs.

IF EXPOSED: Call a POISON CENTER or doctor/physician. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Keep away from heat, sparks, open flames, and hot surfaces. No Smoking. Do not breathe dust/fumes/gas/mist/vapors/spray.

ANTIDOTE
Ingestion: If victim is conscious and alert, give 2-4 cupfuls of milk or water. Never give anything by mouth to an unconscious person. Get medical aid immediately. Induce vomiting by giving one teaspoon of Syrup of Ipecac.


THIS PRODUCT CONTAINS METHYL ALCOHOL (Methanol). Methanol is a poison containing harmful vapors. It may be fatal or cause blindness if swallowed and cannot be made non-poisonous. Refer to the SDS for additional information.

UN 1987
Alcohols N.O.S., (Ethanol, Isopropanol) 3, II

Ethanol CAS No. 64-17-5
Isopropanol CAS No. 67-63-0
Methanol CAS No. 67-58-1


This product meets and exceeds all grade specifications as listed in ACS Reagent Catalog. Refer to the Certificate of Analysis for the actual results for this lot.

84029-ALCO-100-REV03



ReOrder# R54031 32 Commerce Street, Springfield, NJ 07081
t 888.392.8411 www.avantik-us.com

ALCOHOL 95%



**Reagent
95% Alcohol**
190 PROOF REAGENT GRADE

Formula (Denatured 3A)	190 Proof	Percent by Volume
SDA 3A (Ethanol-Methanol)	95 parts by volume	Ethanol 85.95%
IPA	5 parts by volume	Methanol 4.52%
		IPA 5%
		Water 4.52%

1 Gallon | Case = 4 x 1 Gallon

DANGER! Highly flammable liquid and vapor. Harmful if swallowed. Causes skin irritation. Causes serious eye irritation. May cause respiratory irritation. Causes damage to organs.

IF EXPOSED: Call a POISON CENTER or doctor/physician. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Keep away from heat, sparks, open flames, and hot surfaces. No Smoking. Do not breathe dust/fumes/gas/mist/vapors/spray.

ANTIDOTE
Ingestion: If victim is conscious and alert, give 2-4 cupfuls of milk or water. Never give anything by mouth to an unconscious person. Get medical aid immediately. Induce vomiting by giving one teaspoon of Syrup of Ipecac.

THIS PRODUCT CONTAINS METHYL ALCOHOL (Methanol). Methanol is a poison containing harmful vapors. It may be fatal or cause blindness if swallowed and cannot be made non-poisonous. Refer to the SDS for additional information.

UN 1987
Alcohols N.O.S., (Ethanol, Isopropanol) 3, II

Ethanol CAS No. 64-17-5
Isopropanol CAS No. 67-63-0
Methanol CAS No. 67-58-1
Water CAS No. 7732-18-5

Refer to the Certificate of Analysis for the actual results for this lot.

84031-ALCO-95-REV03



ReOrder# R54078 32 Commerce Street, Springfield, NJ 07081
t 888.392.8411 www.avantik-us.com

ALCOHOL 80%



**Reagent
80% Alcohol**
160 PROOF REAGENT GRADE

Formula (Denatured 3A)	160 Proof	Percent by Volume
SDA 3A (Ethanol-Methanol)	95 parts by volume	Ethanol 72.40%
IPA	5 parts by volume	Methanol 5.60%
		IPA 4.00%
		Water 20%

1 Gallon | Case = 4 x 1 Gallon

DANGER! Highly flammable liquid and vapor. Harmful if swallowed. Causes skin irritation. Causes serious eye irritation. May cause respiratory irritation. Causes damage to organs.

IF EXPOSED: Call a POISON CENTER or doctor/physician. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Keep away from heat, sparks, open flames, and hot surfaces. No Smoking. Do not breathe dust/fumes/gas/mist/vapors/spray.

ANTIDOTE
Ingestion: If victim is conscious and alert, give 2-4 cupfuls of milk or water. Never give anything by mouth to an unconscious person. Get medical aid immediately. Induce vomiting by giving one teaspoon of Syrup of Ipecac.


THIS PRODUCT CONTAINS METHYL ALCOHOL (Methanol). Methanol is a poison containing harmful vapors. It may be fatal or cause blindness if swallowed and cannot be made non-poisonous. Refer to the SDS for additional information.

UN 1987
Alcohols N.O.S., (Ethanol, Isopropanol) 3, II

Ethanol CAS No. 64-17-5
Isopropanol CAS No. 67-63-0
Methanol CAS No. 67-58-1
Water CAS No. 7732-18-5


Refer to the Certificate of Analysis for the actual results for this lot.

84078-ALCO-80-REV03



ReOrder# R54033 32 Commerce Street, Springfield, NJ 07081
t 888.392.8411 www.avantik-us.com

ALCOHOL 70%



**Reagent
70% Alcohol**
140 PROOF REAGENT GRADE

Formula (Denatured 3A)	140 Proof	Percent by Volume
SDA 3A (Ethanol-Methanol)	95 parts by volume	Ethanol 63%
IPA	5 parts by volume	Methanol 3.2%
		IPA 3.1%
		Water 30%

1 Gallon | Case = 4 x 1 Gallon

DANGER! Highly flammable liquid and vapor. Harmful if swallowed. Causes skin irritation. Causes serious eye irritation. May cause respiratory irritation. Causes damage to organs.

IF EXPOSED: Call a POISON CENTER or doctor/physician. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Keep away from heat, sparks, open flames, and hot surfaces. No Smoking. Do not breathe dust/fumes/gas/mist/vapors/spray.

ANTIDOTE
Ingestion: If victim is conscious and alert, give 2-4 cupfuls of milk or water. Never give anything by mouth to an unconscious person. Get medical aid immediately. Induce vomiting by giving one teaspoon of Syrup of Ipecac.

THIS PRODUCT CONTAINS METHYL ALCOHOL (Methanol). Methanol is a poison containing harmful vapors. It may be fatal or cause blindness if swallowed and cannot be made non-poisonous. Refer to the SDS for additional information.

UN 1987
Alcohols N.O.S., (Ethanol, Isopropanol) 3, II

Ethanol CAS No. 64-17-5
Isopropanol CAS No. 67-63-0
Methanol CAS No. 67-58-1
Water CAS No. 7732-18-5

Refer to the Certificate of Analysis for the actual results for this lot.

84033-ALCO-70-REV03

Call today and learn how Avantik's color-coded system can eliminate reagent mix-ups and help you...

Achieve TOTAL **WORKFLOW CONTROL**SM



avantik-us.com
800-783-9424



SAVE YOUR SPECIMEN BY ELIMINATING REAGENT MIX-UPS.

Match Before You Process.

Avantik's Matched Tissue Processor Bottle Labels



Call today and learn how Avantik's color-coded system can eliminate reagent mix-ups and help you...

Achieve TOTAL **WORKFLOW** CONTROLSM

avantik-us.com
800-783-9424

