



The embedding center that you choose plays a critical role in the performance of your lab. Throughput, ease of use, ergonomics and performance are all key factors for your lab in order to run efficiently and flawlessly. The Avantik EM8 Embedding Center is an integral part of creating a lab. The modular design allows for easy and customized setup of the hot and cold plate surfaces. The control panel is intuitive and simple.

Features & Benefits

- Ergonomically designed to enhance operator comfort
- Adjustable paraffin flow rate ensures operator efficiency
- Individual control of all temperature settings
- User friendly keypad for easy control of all system operations
- Large cassette storage area accommodates multiple sized tissue baskets
- 5 liter paraffin reservoir ensures continuous workflow
- Large cooling surface holds 60+ cassettes and molds

Product Information

EM8H Dispensing Console

EM8C Cryo Console

Related Accessories

EQ6109 Electric Heated Forceps - Red - 1mm Tip

EQ6110 Electric Heated Forceps - Yellow - 2mm Tip

EQ6111 Electric Heated Forceps - Orange - 4mm Tip

Call us today.

800-783-9424 info@avantik-us.com avantik-us.com





EM8H - Dispensing Console Specifications

 Outside dimensions
 23.25"W x 17"H x 26"D

 Weight
 62 lbs.

 Paraffin reservoir
 5 Liters

 Left & right trays capacity
 1.7 Liters

 Tray measurements
 7"W x 10.5"D x 2.25"H

 Paraffin reservoir temperature range
 50 - 70°C

 Left tray temp
 50 - 70°C

 Right tray temp
 50 - 70°C

 Working area
 50 - 70°C

 Forceps receptacle
 50 - 75°C

 Power supply
 100-240V / 50-60Hz

EM8C - Cryo Console Specifications

Call us today.

800-783-9424 info@avantik-us.com avantik-us.com

