

FLAMMABLE STORAGE CABINET



Our safety storage cabinets serve several critical functions. The greatest of which is to safely contain hazardous chemicals to reduce the risk of fire, protecting both personnel and facilities. Another function of these cabinets is to identify, organize, and segregate dangerous liquids. These cabinets offer compliance with federal OSHA regulations and state and local fire codes. They also improve efficiency by allowing flammable liquids to be stored near the points of use.

Features

- All welded, double wall, 18 gauge construction
- Compliant with federal OSHA regulations and state and local fire codes
- Manual close doors
- Twin 2" flame arrester vents
- 2" leakproof sill to contain leaks
- Lockable flush mounted handle (gripper pad) with 2 keys
- Adjustable leveling feet and grounding connector
- Galvanized steel shelves adjust on 2-1/2" centers

Benefits

- Identify, organize, and segregate dangerous liquids
- Reduce risk of fire
- Protect personnel and facilities

Product Information

Yellow Flammable Liquid Safety Cabinets

Item	Description	Dimensions	Weight
SY5095	4 gallon	17" W x 17" D x 22"H	63 lbs
SY6068	12 gallon	23.25" W x 18" D x 35.125"H	150 lbs
SY6019	15 gallon	23.25" W x 18" D x 44" H	167 lbs
SY5093	24 gallon	23" W x 16" D x 65" H	194 lbs
SY6025	44 gallon	43" W x 18" D x 65" H	345 lbs
SY6016	45 gallon	43" W x 18" D x 65" H	358 lbs
SY6024	55 gallon	34" W x 34" D x 65" H	350 lbs
SY5090	60 gallon	34" W x 34" D x 65" H	430 lbs
SY6067	90 gallon	43" W x 34" D x 65" H	461 lbs

White Flammable Liquid Safety Cabinets

Item	Description	Dimensions	Weight
SY5030	4 gallon	17" W x 17" D x 22"H	63 lbs
SY5032-B	12 gallon	23.25" W x 18" D x 35.125"H	150 lbs
SY5032-C	22 gallon	35" W x 22" D x 35" H	170 lbs
SY5032-D	30 gallon	43" W x 18" D x 44" H	240 lbs
SY5032-E	45 gallon	43" W x 18" D x 65" H	358 lbs
SY5032-F	60 gallon	34" W x 34" D x 65" H	430 lbs

Piggyback Flammable Liquid Safety Cabinets

Item	Description	Dimensions	Weight
SY6203-A	12 gallon	23.25" W x 18" D x 35.125"H	129 lbs
SY6203-B	15 gallon	45" W x 22" D x 29" H	169 lbs

Call us today.

800-783-9424
info@avantik-us.com
avantik-us.com

