


THE AVANTIK **LP16**  
 LABEL PRINTING SYSTEM



Thousands of specimen identification errors occur annually which inevitably leads to adverse patient events. A large percentage of these errors are due to illegible slides. Combined with Avantik's optimized Mohs Software Solution, the LP16 printer brings the confidence and reliability of accurate specimen identification to your laboratory.

### Key Benefits

- Customizable Mohs specific software allows for fast data entry and tech efficiency.
- Chemical resistant labels which means that the slide can be stained with no risk of losing patient data.
- Daily Log File creates a record of every case printed and creates CLIA log file eliminating handwritten CLIA logs.

### Product Information

**LP16** LP16 Printer, Printer Ribbon, Labels, USB Cable, PC Touchscreen, Mohs Software, Monitor  
 \*VESA Mountable

### Printer Specifications

**Width**.....7" or 177.8 mm  
**Height**.....8" or 203.2 mm  
**Depth**.....10" or 254 mm  
**Weight**.....4.6 lbs. or 2.09 kg.

### Screen Specifications

**Screen Dimensions** .....15.6" or 396.24 mm

Call us today.

**800-783-9424**  
**info@avantik-us.com**  
**avantik-us.com**

