

SP15

SLIDE PRINTER



Thousands of specimen identification errors occur annually, which inevitably leads to adverse patient events. A large percentage of these errors are due to improperly labeled slides. The small footprint of the new Avantik SP15 brings the convenience of on-demand slide printing directly to your cutting station. Now, with the Avantik SP15 Slide Printer, you can eliminate errors by allowing the technicians to reliably print their slides at the time the specimens are cut.

Features

- Small footprint fits easily next to the microtome
- Chemical resistant media ensures print reliability
- Optional unload carrier
- Flexible data inputs allow easy LIS integration
- Comes in two different models

Benefits

- Eliminates the need for labels or handwriting
- Media ensures print reliability
- Easy integration with the LIS system
- Slides protected from dust & the environment
- Enhanced productivity in the lab

Product Information

SP15 SP15 Slide Printer

SP15-U SP15 Slide Printer with Unload Station

Accessories

SL6342 SP15 Slide Printer Ribbon

SL6332 SureBond Charged Microscope Slides

SL6333 SureSlide Standard Microscope Slides

Call us today.

800-783-9424
info@avantik-us.com
avantik-us.com



NEW AVANTIK
SP15
 SLIDE PRINTER



SP15 Specifications

Width	6.0", 152mm
Height.....	8.9", 226mm
Depth.....	10.4", 264mm
Weight.....	4.3kg, 9.5 lbs.
Power Requirements	Universal 100-240 VAC 47-63 Hz 1.2A Maximum
Color Bone White	
Slide Type.....	3" x 1" or 25mm x 75mm x 1mm
Slide Capacity	150 pieces (input stack) 10 pieces (output tray)
Media Type.....	Thermal transfer ribbon, 1" wide roll
Media Capacity.....	1 roll. Approx. 12,000 Slides
Media Color.....	Black
Throughput.....	8 slides/minute typical

Call us today.

800-783-9424
info@avantik-us.com
avantik-us.com

