


TrimSafe SERIES



The Avantik TrimSafe Series further reduces the risk of injury from blade malfunction and user error. The TrimSafe handle is large for improved ergonomics and it is easy to clean. The TrimSafe blades have a precision edge for long lasting sharpness. The push lever removes blades quickly for convenient disposal. The series also comes with a TrimSafe dispenser which makes it easy and safe to load new blades.



TrimSafe HANDLE

- Quicker and safer loading of trim blades onto handle
- Larger molded handle for improved ergonomics
- Push lever quickly clamps to secure and release TrimSafe blades

130mm Handle - Item# GR2164-A
260mm Handle - Item# GR2164-B



TrimSafe BLADES

- Precision sharpened edge for ultimate sharpness and durability
- Back edge of blade rounded to eliminate tissue snags and to improve user's safety

130mm Blades - Item# GR2163-A
260mm Blades - Item# GR2163-B



TrimSafe DISPENSER

- Load & unload trim blades without the need to touch the blades
- Pre-loaded with premium TrimSafe Blades for ease and compatibility
- Built-in sharps container for easy and safe disposal of used blades

130mm Dispenser - Item# GR2173-A
260mm Dispenser - Item# GR2173-B

Call us today.

800-783-9424
info@avantik-us.com
avantik-us.com

