

AVANTIK ULTRAFIN®

A SOLUTION FOR EVERY LAB



Every laboratory has different requirements for their paraffin and **Avantik UltraFin®** exceeds them all. One brand, three unique solutions for every lab. These three paraffin solutions offer great infiltration, easy ribboning, and are in stock, ready to ship.



Avantik UltraFin® Paraffin

This tissue embedding medium with rapid penetration is suggested for large tissue or tissue that may present infiltration issues. Promotes excellent ribbon continuity to 2um; melts quickly at **56-57°C**. The double-filtered paraffin contains plastic polymers and (0.8%) dimethyl sulphoxide (DMSO) for faster tissue penetration. Compare to Paraplast Plus™.

Item	Size	Qty
EM0125-A	2.2 lb. Bag	8/Case



Avantik UltraFin® IE Paraffin

The ultimate infiltration medium is recommended for general tissue embedding; no more icing blocks are needed. This paraffin reduces section wrinkles and gives difficult to section tissue more elasticity. The medium supports improved tissue infiltration and quality sections with minimal compression. Melts at **56-58°C**. Compare to Paraplast™.

Item	Size	Qty
EM0125-B	2.2 lb. Bag (1kg)	8/Case



Avantik UltraFin® EX Paraffin

Avantik UltraFin® EX is recommended for lower temperature infiltration, which aids in the prevention of tissue distortion. This paraffin is ideal for working with heat-sensitive tissue. The medium's lower viscosity allows complete infiltration of dense tissue. Melts at **50-54°C**. Compare to Paraplast X-tra™.

Item	Size	Qty
EM0125-C	2.2 lb. Bag (1kg)	8/Case

Peace of Mind,
Guaranteed Every Day.
Avantik, **A Higher Standard.**



Avantik UltraFin[®] **A Solution for Every Lab**

One brand, three unique solutions. Avantik UltraFin[®] Paraffin solutions offer great infiltration, easy ribboning, and are in stock, ready to ship.



SCAN NOW

to learn more about
Avantik UltraFin[®]

Contact Us Today!
avantik-us.com
800.783.9424

Peace of Mind,
Guaranteed Every Day.
Avantik, **A Higher Standard.**