

THE AVANTIK EM8 EMBEDDING CENTER

AN INTEGRAL PART OF YOUR LAB EFFICIENCY



The embedding center that you choose plays a critical role in the performance of your lab. The modular design of the Avantik EM8 Embedding Center allows for a customized setup to help elevate efficiency in your lab.

User Benefits

- Ergonomically designed to enhance operator comfort
- Cryo console operates on its own, which allows use on either side of the embedding console
- User friendly keypad for easy control of all system operations
- Adjustable paraffin flow rate ensures operator efficiency
- Individual control of all temperature settings

Product Details

- Large cassette storage area accommodates multiple sized tissue baskets
- 5 liter paraffin reservoir ensures continuous workflow
- Large cooling surface holds 60+ cassettes and molds

Scan Now to View
Product Video



Peace of Mind.
Guaranteed Every Day
Avantik, **A Higher Standard.**



THE AVANTIK

EM8 EMBEDDING CENTER

AN INTEGRAL PART OF YOUR LAB EFFICIENCY

EM8H - Dispensing Console

Specifications

Outside dimensions	23.25"W x 17"H x 26"D
Weight	62 lbs.
Paraffin reservoir	5 Liters
Left & right trays capacity	1.7 Liters
Tray measurements	7"W x 10.5"D x 2.25"H
Paraffin reservoir temperature range	50 - 70°C
Left tray temp.	50 - 70°C
Right tray temp.	50 - 70°C
Working area	50 - 70°C
Forceps receptacle	50 - 75°C
Power supply	100-240V / 50-60Hz

EM8C - Cryo Console

Specifications

Outside dimensions	13"W x 16"H x 23.75"D
Cooling surface.....	12" x 14.75"
Weight	55 lbs.
Cooling unit.....	0 to -15°C

Related Accessories

- EQ6109** Electric Heated Forceps - Red - 1mm Tip
- EQ6110** Electric Heated Forceps - Yellow - 2mm Tip
- EQ6111** Electric Heated Forceps - Orange - 4mm Tip

Contact Us Today!

avantik-us.com
800.783.9424

