

ESPO SLIDE PRINTER



The ESPO Slide Printer is easy to use, provides efficiency for your lab, and has the flexibility to print on slides containing mounted tissue.

It revolutionizes the laboratory by providing on-demand printing, creating a more efficient workflow and minimizing errors.

Specifications

Capacity	Two 100-slide magazines, capable of printing 2 different types or colors of slides
Control Slide	Single-slide front loading for control slides with tissue on slide
Speed	Print speed is 3 seconds per slide / 20 per minute
Ease of Use	Fully integrated color touch-screen display and simple software is easy to use
Small Footprint	Ideal for on-demand printing

Features

- Fully integrated color touchscreen display and simple software is easy to use
- Two 100-slide magazines allow for continuous printing
- Single slide front loading enables refill-free printing, even when tissue is on the slide

Product	Order Number
ESPO Slide Printer	EQ6142
50 Slide Stacker	EQ6152
Ribbon	SL6508
Slide Magazine	SL6509

Call us to order!

avantik-us.com
800.783.9424

