

ESPO SLIDE PRINTER

A HIGHER STANDARD™ FOR SLIDE LABELING

The **ESPO Slide Printer** is an easy-to-use, compact-designed, thermal-transfer printer designed to provide efficiency for your lab. The ESPO has the flexibility to print on slides with or without mounted tissue. Say goodbye to handwritten labels and identification errors with the confidence provided by ESPO.

Product Features

- 1 Display Screen
- 2 Two 100-Slide Cartridges
- 3 Small Footprint




User Benefits



Simplified Operation
Fully integrated, icon based, touch-screen interface is easy to use.



Slide Capacity
Two 100-slide cartridges allow for continuous printing; single slide front loading option for refill-free printing.




Easy Connectivity
You can connect the ESPO printer to a desktop or laptop computer, or to a LIS via WiFi, LAN or USB.



Compact Design
The ESPO takes up minimal space alone or within the microtome benchtop space.

ESPO Printer Service Difference

With the **ESPO Slide Printer**, Avantik® offers a full range of services, including **on-site installation, training and support.**



Peace of Mind,
Guaranteed Every Day.
Avantik, **A Higher Standard™**



ESPO Printer Specifications

Slide Capacity	Load up to 200 slides via two, 100-slide magazines
Printing Speed	Print time is 3 seconds per slide, 20 slides per minute
Printing Type	Thermal transfer ribbon, 12 dot/mm
Dimensions (D x W x H)	9.4" x 4.7" x 9.8"
Small Footprint	15.4 lbs.

Ordering Information

Inquire about our ESPO Printer service offerings

EQ6142	ESPO Slide Printer
EQ6152	50 Slide Stacker
SL6508	Ribbon
SL6509	Slide Magazine



SCAN NOW
to Request a
Virtual Demo

Contact Us Today!
avantik-us.com
800.783.9424

MAR-0000189 Rev. 1

 **AVANTIK**[®]
A HIGHER STANDARD[™]

