

Quality. Consistency. Reliability. **The New Avantik LP16 Label Printing System.**

NEW!
LP16 LABEL PRINTER

THE AVANTIK

LP16

LABEL PRINTING SYSTEM

Your Lab...Quick, Efficient, Easy

The New Avantik LP¹⁶ Label Printer

Thousands of specimen identification errors occur annually, which inevitably leads to adverse patient events. A large percentage of these errors are due to illegible slides. Combined with Avantik's optimized Mohs Software Solution, the LP¹⁶ brings the confidence and reliability of accurate specimen identification to your laboratory.

Benefits of the Avantik LP¹⁶ Label Printing System

- *Eliminate slide identification errors*
- *Small footprint fits anywhere in your lab*
- *Chemical resistant labels to ensure high quality and permanency*
- *Cost-effective solution*
- *Optimized software solution for Mohs*
- *Fast and reliable slide labeling system for improved lab efficiency*



THE MOHS SOLUTION TO ACCURATE AND DEPENDABLE SPECIMEN IDENTIFICATION

Call Now To Order!

www.avantik-us.com

888.392.8411

 **AVANTIK**

THE AVANTIK LP16

LABEL PRINTING SYSTEM

NEW!
LP16 LABEL PRINTER



SPECIFICATIONS

- Width.....7", 177.8mm
- Height.....8", 203.2mm
- Depth.....10", 254mm
- Weight.....2.09kg. 4.6 lbs.
- Power RequirementsUniversal 100-240 VAC
47-63 Hz 1.2A Maximum
- Color.....Charcoal
- Print Resolution.....300 DPI
- Maximum Print Width.....4"

ORDER INFORMATION

- LP16 Label Printer**.....LP16
Includes LP16 Printer
- LP16 Label Printing System**.....LP16-K
Includes LP16 Printer, Mohs Software, Printer Ribbon, Labels, USB Cable

THE MOHS SOLUTION TO ACCURATE AND DEPENDABLE SPECIMEN IDENTIFICATION

Call Now To Order!
www.avantik-us.com
888.392.8411

